DENTAL MACHINE

CAM – T15

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INSTRUCTION FOR USE INDICATIONS FOR USE, DESCRIPTION AND COMPOSITION OF THE MEDICAL DEVICE

The Medical Device CAM T15 Disc in Titanium alloy for dentistry is codified as follows: cd-xx (xx may vary from 01 to 11 based on the thickness of the blank, the instruction for use remains the same for all the variants). The disc is composed by a dental alloy made of Titanium conformed to ISO 5832-3_2016, intended to be milled by CNC machines, for the manufacturing of dental prosthetic elements.

CAM-T15 can be used to manufacture components with thin sections exposed to high stresses like esthetic coated crowns, extended bridges or bridges with small sections, bars, anchors, abutments.

CAM-T15 is intended for professional use only, specifically for dentists and dental technicians.

The alloy composition is as follows:

Chemical Element	Composition % [m/m]
Aluminium	5.5 – 6.75
Vanadium	3.5 – 4.5
Iron	0.3 max
Oxygen	0.2 max
Carbon	0.08 max
Nitrogen	0.05 max
Hydrogen	0.015 max
Titanium	Remaining part

The alloy is completely conformed to ISO 5832-3:2016's regulations "Implants for surgery — Metallic Materials, Part 3: Wrought titanium 6-aluminium 4-vanadium alloy".

MODES OF ACTION AND PERFORMANCE

The material has to be milled by a CAD-CAM-software-based CNC machine in order to manufacture several elements in the dental sector.

Once the piece has been designed with the milling machine software, the disc has to be placed in the disc holder and fastened to avoid any shift during the milling process.

The elements can be modeled and/or finished using proper tools for the manufacturing of titanium (ATTENTION: sides and edges must not be sharp).

It is advisable to direct all the tools to the surface to one direction only, in order to avoid the overlapping of the material. Moreover, you should not overpass the number of rotations recommended by the tool's manufacturer.

Before proceeding with any sort of esthetic coating, you have to clean all the components, the procedure is left to the user's craftsmanship.

PACKAGING FORMAT

The pack contains one single disc in Titanium CAM – T15 alloy.

CONDITIONS OF USE AND HANDLING

Titanium alloys are intrinsically high resistant to corrosion thanks to the Titanium dioxide layer that becomes visible immediately after the manufacturing. Although, you should keep the medical device in its own pack until use, protected from the sunlight and moisture, because damage may occur on the label information as well as in the material (in very rare cases).

PRECAUTIONS

- Conform to the instruction of use
- Respect the intended use provided and described in the Medical Device description
- A hand-made finishing of the component may be necessary after the milling process
- The element can be coated with resins and composites (in this case, you must conform to the instruction provided by the material's manufacturer)
- The Medical Device has to be handled by qualified staff (dental technicians and dentists)
- The final element has to be cleaned before use
- Do not melt the disc residuals
- Do not expose the material to temperatures higher than 800°C, variations in the microstructure may occur
- The metal steams and powder are harmful for the health in case of inhalation; therefore, a suction system or suitable protection masks must always be employed during the manufacturing process.
- The risk of fire subsists in the manufacturing of Titanium and Titanium-based alloys.
- The CNC milling machine has to be used with constant use of lubricant/cooling liquid.
- Do not ever use Hydrofluoric acid.

COLLATERAL EFFECTS AND CONTRINDICATIONS

The Device is a dental alloy and does not present any contraindications if it is employed for the intended use provided and described in the instruction.
The Device must NOT be employed to manufacture any elements that are not intended for the dental sector because the manufacturer has NOT taken into consideration other medical sectors.

STORAGE CONDITIONS

The device must be kept away from the direct sunlight and protected from moisture.

LIMITATIONS OF USE

The Medical Device is intended for the use described in the paragraph "INDICATIONS FOR USE, DESCRIPTION AND COMPOSITION OF THE MEDICAL DEVICE". Any other use is forbidden, especially the manufacture of elements that are not intended for the dental sector.

DISPOSAL

The material residuals from the manufacturing and any elements manufactured using the Medical Device must be disposed of in conformity with the national legislation.

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